## SERVATOR® B

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DATA DI REVISIONE DEL TESTO: Marzo 2023

## SERVATOR® B

STERILE AND PYROGEN-FREE SOLUTION FOR ORGAN FLUSHING AND PRESERVATION  
NOT INJECTABLE

**COMPOSITION:** 1000 ml of solution contain:

Poly (o-2-hydroxyethyl) starch 0.40-0.50 MS <sup>1</sup>	(50.0 g/l)	
Lactobionic acid (as lactone)	(35.83 g/l)	105 mmol/l
Potassium hydroxide	(5.61 g/l)	100 mmol/l
Sodium hydroxide	(1.08 g/l)	27 mmol/l
Adenosine	(1.34 g/l)	5 mmol/l
Allopurinol	(0.136 g/l)	1 mmol/l
Potassium dihydrogen phosphate	(3.4 g/l)	25 mmol/l
Magnesium sulphate heptahydrate	(1.23 g/l)	5 mmol/l
Raffinose pentahydrate	(17.83 g/l)	30 mmol/l
Glutathione	(0.922 g/l)	3 mmol/l
Potassium hydroxide/Phosphoric acid	to pH adjustment q.s. to pH 7.4	
Water for injection	q.s. to 1000 ml	

<sup>1</sup>MS = moles of hydroxyethyl groups per moles of anhydroglucose units

**PHYSICAL PROPERTIES:**

The solution has a calculated osmolarity of 320 mOsm/l, a sodium ions final concentration of 29 mEq/l, a potassium ions final concentration of 125 mEq/l.  
pH of 7.1 - 7.5 at 20°C.

**TYPE OF DEVICE:**

Sterile and pyrogen-free solution for organ preservation. Clear, colorless, or slightly yellow solution.

**INDICATIONS:**

Flushing and hypothermic preservation of kidney, liver, and pancreas.

At a temperature of 2°C-6°C (36°F-43°F), Servator® B solution is used to flush the isolated organ immediately before or immediately after removal from a dead donor or immediately after removal from a living donor.

The solution remains in the organ vessels during hypothermic storage and transportation.

Servator® B solution is meant for the cold storage of the organ and is not indicated for hypothermic preservation with continuous mechanical perfusion.

The use of the solution at the recommended temperatures allows an effective cooling of the organ.

The device should be used only by medical staff properly trained according to established operating protocols.

**CONTRAINDICATIONS:**

Hypersensitivity to one of its components, in particular to Allopurinol and Poly (o-2-hydroxyethyl) starch 0.40-0.50, or to any of the additives (penicillin, insulin and dexamethasone).

**SPECIAL WARNINGS AND PRECAUTIONS FOR USE:**

NOT FOR DIRECT INJECTION OR INTRAVENOUS INFUSION.

Systemic administration may cause side effects that can endanger the patient's life.

**CAUTION:**

As Servator® B solution contains a high concentration of potassium, flush the donor organ free of residues with a physiological solution before reperfusion during transplantation, to avoid severe cardiovascular complications such as cardiac arrest due to hyperkalaemia or bradyarrhythmia.

The incomplete elimination of residues of the preservation solution may result in hypersensitivity to any of its components or additives (penicillin, insulin, and dexamethasone).

**INSTRUCTIONS FOR USE:**

Cool Servator® B solution to a temperature of 2°C-6°C (36°F-43°F). Remove the overwrap before use. Do not use if the bag under pressure leaks or if visible particles, precipitate or contamination are detected upon visual inspection.

Do not filter before use (the solution undergoes 0.2 micron filtration during production process).

Immediately prior to use and under aseptic conditions add to each liter of Servator® B solution the following drugs:

200.000 units of penicillin G;

40 units of rapid insulin;

16 mg dexamethasone.

The glutathione present in the solution is subject to oxidation during preservation of the organs. Consequently, 0.922g/l (3 mmol/l) of freshly prepared solution of reduced glutathione may be added immediately before use.

Remove the protective cap from the outlet fitting and insert the administration set spike with a rotary motion. Open the clamp of the administration set. Keep the administration set vertically above the solution bag and compress the bag to fill the administration set with fluid. Close the clamp.

Prior to connection to the organ, suspend the SALF Servator B bag at such a height as to allow a constant flow rate of at least 30 ml/min of solution during flushing.

Open the clamp and start flushing. Continue until the organ is uniformly pale and the output liquid is relatively clear.

Minimum recommended volumes	
In situ aortic flush	Adults 2 – 4 l
	Children 50 ml/kg
Ex vivo infusion Liver (via portal vein and biliary tree)	Adults 1200 ml
	Children 50 ml/kg
Pancreas or kidney	Adults 300 – 500 ml
	Children 150 – 250 ml

Add more solution in the storage container until the organ is fully immersed. Seal the container aseptically and place it in another container well insulated and suitable for organ transportation.

Place ice on the outside of the container but not inside, as the ice could enter into direct contact with the organ.

During transplantation and prior to anastomosis, the donor organ should be flushed thoroughly. To reduce residues of solution in the liver flush immediately before the anastomosis with one liter of lactated Ringer’s solution through the portal vein.

**ISCHEMIC TIME:**

The recommended time for each organ is the following:

<i>Cold ischemia time</i>	
Liver	Not more than 17 hours
Kidney	Not more than 23 hours
Pancreas	Not more than 21 hours

<i>Warm ischemia time</i>	
---------------------------	--

Liver	Not more than 2.5 hours
Kidney	Not more than 2.5 hours
Pancreas	Not more than 2.5 hours

**INCOMPATIBILITY:**

No known interactions when used as indicated.

Servator® B solution should not be combined with medicinal products except the drugs referred to in paragraph "Instructions for use".

**SIDE EFFECTS:**

If the cold storage solution is not completely eliminated from the organ after flushing with a physiological solution and prior to transplantation, the recipient may undergo cardiovascular complications, like bradyarrhythmia.

Clinical problems have been reported in relation to the conservation of the liver including impairment of liver function, liver failure, death of the patient.

Histologically, ischemic lesions have been observed in association with or even in absence of cases of mild rejection.

**WARNINGS AND PRECAUTIONS IN RELATION TO THE SHELF LIFE OF THE DEVICE:**

Check the expiration date on the container. Shelf life refers to the product in its intact packaging and correctly stored.

CAUTION: Do not use after the expiration date. Do not use the solution even before expiration if visible particles, precipitate or contamination are detected.

Use the cold storage solution immediately after addition of the additives.

**SPECIAL PRECAUTIONS FOR STORAGE:**

Sterile and disposable device, for one single and continuous administration. Discard any residue to avoid risk of contamination due to loss of sterility.

Store at a temperature of 2°C-25°C (36°F-77°F). Avoid exposure to heat. Do not freeze.

Do not use if the solution is frozen. Store in its original container. Do not remove the overwrap until immediately prior to use.

KEEP OUT OF REACH OF CHILDREN.

**TYPE OF PACKAGING:**

PVC-free bag containing 1000 ml of solution.

Box of 10 bags.

PVC-free bag containing 2000 ml of solution.

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**SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP):**

The document is available at the following URL: <https://cloud.salfspa.it:8443/index.php/s/SSCP.ServatorH>

**SERIOUS INCIDENT REPORTING:**

The user and/or patient are recommended to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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S.A.L.F. S.p.A. LABORATOIRE PHARMACOLOGIQUE Via Marconi, 2 24069 Cenate Sotto, Bergamo - Italie