SERVATOR® C

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DATA DI REVISIONE DEL TESTO: Marzo 2023

SERVATOR® C

STERILE AND PYROGEN-FREE SOLUTION FOR ORGAN FLUSHING AND PRESERVATION
NOT INJECTABLE

COMPOSITION: 1000 ml of solution contain:

Sodium hydroxide	4.0 g	100 mmol/l
Potassium chloride	1.12 g	15 mmol/l
Calcium chloride dihydrate	0.037 g	0.25 mmol/l
Magnesium chloride Hexahydrate	2.645 g	13 mmol/l
L-Histidine	4.655 g	30 mmol/l
Lactobionic acid	28.66 g	80 mmol/l
Mannitol	10.931 g	60 mmol/l
Reduced glutathione	0.922 g	3 mmol/l
Glutamic acid	2.943 g	20 mmol/l
Water for injections q.s. to 1000 ml		

PHYSICAL PROPERTIES:

The solution is slightly hypertonic, with a calculated osmolarity of about 242-368 mOsmol/l, low viscosity and high buffer property.

Slightly alkaline, pH 7.1 - 7.5 at 20°C.

TYPE OF DEVICE:

Sterile and pyrogen-free solution for organ preservation. Clear, yellow solution.

INDICATIONS:

Flushing and hypothermic preservation of thoracic organs (heart and lung) as well as abdominal organs (kidney, liver, and pancreas) during transplantation.

At a temperature of 2°C-8°C (36°F-46°F), Servator® C solution is used to flush the isolated organ immediately before or immediately after removal from a dead donor or immediately after removal from a living donor. The solution remains in the organ vessels during hypothermic storage and transportation. Servator® C solution is meant for the cold storage of the organ and is not indicated for hypothermic preservation with continuous mechanical perfusion.

The use of the solution at the recommended temperature allows the proper cooling of the organ.

Do not filter before use (the solution undergoes 0.2 micron filtration during production process).

The device should be used only by medical staff properly trained according to established operating protocols.

CONTRAINDICATIONS:

Hypersensitivity to one of its components.

SPECIAL WARNINGS AND PRECAUTIONS FOR USE

NOT FOR DIRECT INJECTION OR INTRAVENOUS INFUSION.

Servator® C solution is to be used only for organ flushing and preservation during transplantation. Its safety and effectiveness in preserving organs from pediatric donors or to be transplanted in pediatric patients have not been studied.

- Do NOT use in vivo
- Do NOT use for direct injection
- Do NOT use for intravenous infusion
- Do NOT use for continuous perfusion
- The solution is sterile and disposable.

- During transplantation and prior to reperfusion, flush the donor organ accurately with a physiological solution to avoid severe cardiovascular complications such as cardiac arrest due to hyperkalaemia or bradyarrhythmia. The incomplete elimination of residues of the preservation solution may result in hypersensitivity to any of its components.
- Do NOT use the solution if frozen
- Do NOT administer at systemic level
- DISCARD the bag if visible particles, precipitate, or contamination are present in the solution
- Do NOT use if the bag under pressure shows leakage
- Do NOT dilute nor add medications or additives
- Do NOT ingest
- DISCARD the solution after the expiration date
- Do NOT filter.

INSTRUCTIONS FOR USE:

- Take the Servator® C solution and cool it to a temperature between 2°C-8°C (36°F-46°F) before use.
- Use only if the bag and the plastic overwrap are intact.
- Remove the plastic overwrap and control the appearance of the solution: it should not be cloudy nor contain precipitates or contamination.
- Do not use the solution if frozen.
- Any yellow coloration and/or slight odor of sulphur do not have effects on the quality or efficacy of the solution.
- Disposable solution to be used immediately after opening.
- The solution does not require filtration.
- Not for use with continuous mechanical perfusion.

Place the solution bag inside a sphygmomanometer of suitable dimensions. Inflate the cuff so as to exercise enough pressure on the fluid. Prior to flushing hang the solution bag at such a height as to allow a constant flow.

Open the clamp and start flushing. Continue until the organ is uniformly pale and the output liquid is relatively clear.

The flushing method depends on the centre and on the possible contemporary removal of more organs. Flushing is often performed in two steps: while the organ in situ and once the organ has been removed.

The organ is flushed by means of a cannula inserted in an artery. For a proper rinsing enough pressure should be secured, to allow a constant flow of the solution. In liver transplants, the biliary tree is usually flushed after removal, before being placed in the container for storage and transportation.

The volume of the flushing solution depends on the organ.

<i>Minimum recommended volumes</i>	
Heart	Adults 1 – 2 liters
	Children 30 ml/kg
Liver	6 – 8 liters
Lung	4 – 6 liters
Kidney	4 – 5 liters
Pancreas	4 liters
Multiorgan	According to the single organs involved

Continue flushing until the organ is uniformly pale and the output liquid is relatively clear.

Add more solution in the container until the organ is fully immersed. Seal the container aseptically and place it in another container well insulated and suitable for organ transportation. Place ice on the outside of the storage container but not inside, as the ice could enter into direct contact with the organ.

INCOMPATIBILITY:

No known interactions when used as indicated. Servator® C solution should not be combined with medicinal products or additives.

SIDE EFFECTS:

If the cold storage solution is not completely eliminated from the organ after flushing with a physiological solution and prior to transplantation, the recipient may undergo cardiovascular complications, like cardiac arrhythmia or hypotension. The incomplete elimination of residues of the preservation solution may result in hypersensitivity to any of its components.

WARNINGS AND PRECAUTIONS IN RELATION TO THE SHELF LIFE OF THE DEVICE:

Check the expiration date on the container. Shelf life refers to the product in its intact packaging and correctly stored.

CAUTION: Do not use after the expiration date. Do not use the solution even before expiration if turbidity, visible particles, precipitate, or contamination are detected.

The product should be used within 24 hours after opening.

SPECIAL PRECAUTIONS FOR STORAGE:

Sterile and disposable device, for one single and continuous administration. Discard any residue to avoid risk of contamination due to loss of sterility. Avoid exposure to heat. Do not freeze. Do not use if the solution is frozen.

KEEP OUT OF REACH OF CHILDREN.

Store in its original container. Do not remove the overwrap until immediately prior to use.

TYPE OF PACKAGING:

PVC-free bag containing 1000 ml of solution. Box of 10 bags.

SPECIAL PRECAUTIONS FOR DISPOSAL OF THE DEVICE:

Unused residues of the solution and waste materials should be disposed of in conformity with the local regulations.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP): The document is at the following URL: <https://cloud.salfspa.it:8443/index.php/s/SSCP.ServatorH>

SERIOUS INCIDENT REPORTING: the user and/or patient are recommended to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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S.A.L.F. S.p.A. LABORATOIRE PHARMARCOLOGIQUE Via Marconi, 2 24069 Cenate Sotto, Bergamo - Italie