SERVATOR® H

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DATA REVISIONE: Novembre 2023

SERVATOR® H

STERILE PYROGEN-FREE SOLUTION FOR ORGAN PRESERVATION
NOT INJECTABLE

CE 0373

HTKE0110GPO
HTKE0120GPO
HTKE0150GPO
HTKV0101GOO

COMPOSITION: 1000 ml solution contain:

sodium chloride	0.8766 g	15.0 mMol
potassium chloride	0.6710 g	9.0 mMol
magnesium chloride · 6 H ₂ O	0.8132 g	4.0 mMol
histidine hydrochloride · H ₂ O	3.7733 g	18.0 mMol
histidine	27.9289 g	180.0 mMol
tryptophan	0.4085 g	2.0 mMol
mannitol	5.4651 g	30.0 mMol
calcium chloride · 2 H ₂ O	0.0022 g	0.015 mMol
α-Ketoglutaric acid	0.1461 g	1.0 mMol
Potassium hydroxide	to pH adjustment	
Water for injection	q.s.	

PHYSICAL PROPERTIES: *pH:* 7.02 - 7.20 at 25°C; 7.40 - 7.45 at 4°C; *Osmolality:* 310 mOsm/Kg

INDICATIONS: The device is intended for preservation of organs for transplant with venous or arterial segments (heart, kidney, liver, pancreas). The solution should be used only by medical staff adequately trained according to established operating protocols.

TYPE OF DEVICE AND MODE OF ACTION: Electrolyte and amino acid aqueous solution for organ preservation. Surface cooling of organs (heart, kidney, liver, pancreas). Protection of vessel grafts.

The solution should be cooled before use according to the specific indications listed below for each type of organ.

The cold solution can be used for flushing isolated organs immediately after their removal from the donor according to the protocols thereto. The solution is then left in the organ vasculature during hypothermic preservation. Servator® H Solution should be used for cold storage of organs and not for a continuous machine perfusion.

CONTRAINDICATIONS: None known at present.

PRECAUTIONS AND WARNINGS: Servator® H solution is not intended for intravenous or intra-arterial administration, but only for selective infusion of the arrested heart or the kidney; for surface cooling and preservation of donor organs during transport from donor to recipient, and for preservation of vessel grafts. The solution must not be used for systemic infusion, direct injection or intravenous infusion. In addition, it is not suitable for reintegration of the hematic volume, electrolytes or aminoacids.

INTERACTION WITH DRUGS AND INCOMPATIBILITIES: Interactions with drugs are unknown. The drugs often administered perioperatively are glycosides, diuretics, nitro derivatives, antihypertensives, beta blockers and calcium antagonists. The Servator® H solution must not be mixed with other drugs.

INCOMPATIBILITIES: None known at present.

HEART

RECOMMENDED SOLUTION TEMPERATURE: 5-8°C

For a correct infusion of isolated organs, the following instructions should be followed:

INFUSION VOLUME: 1 ml of solution per minute and per gram (1 ml/min/g) of estimated heart weight (in adults the heart weight is about 0.5% of body weight; in infants it is approx. 0.6% of body weight).

INFUSION PRESSURE (THE SAME PRESSURE AS IN THE AORTIC ROOT): in adults, place the container 50-70 cm above the heart level, equivalent to 40-50 mm Hg. In infants and young children, place the container 40-50 cm above the heart level, equivalent to 30-40 mm Hg.

INFUSION TIME: Under these dosage and pressure conditions, infusion time should be 6-8 minutes. To ensure homogenous equilibration of the myocardium, this time should never be reduced.

INFUSION TECHNIQUE: Monitor carefully infusion time and height of the container. Alternatively, use a infusion pump and monitor time and pressure in the aortic root.

For a correct preservation of organs for transplantation, the following instructions should be followed:

TRANSPLANTATION: A heart intended for transplantation must be stored in the Servator® H solution (2°C - 4°C) so as to guarantee proper protection up to the time of implantation into the recipient.

KIDNEY

RECOMMENDED SOLUTION TEMPERATURE: 5-8°C

For a correct infusion of isolated organs, the following instructions should be followed:

INFUSION VOLUME: 1.5 ml of solution per minute and per gram (1.5 ml/min/g) of kidney's estimated weight. The normal weight of the kidney in adults is about 150 grams.

INFUSION PRESSURE (renal artery): place the solution container at 120-140 cm above the kidney level, equivalent to 90-110 mm Hg, measured at the tip of the infusion catheter in the renal artery.

INFUSION TIME: Under these dosage and pressure conditions, the infusion time will be 8-10 minutes. This time is necessary to ensure homogenous equilibration of the extracellular space of the kidney (including the interstitial tissue and the renal tubule system). Under no circumstances this infusion time should be reduced.

ADDITIONAL MEASURES: To obtain an optimal protective effect of the Servator® H solution on the kidney, it is important to induce in the patient a fast and effective diuresis before starting infusion adopting pharmacological measures and/or hydration.

INFUSION TECHNIQUE: During hydrostatic infusion, time and height of the container above the kidney should be carefully monitored. Alternatively operate with a infusion pump and monitor time and pressure at the tip of the infusion catheter.

For a correct preservation of organs for transplantation, the following instructions should be followed:

TRANSPLANTATION: Kidneys intended for transplantation should be preserved in the Servator® H solution at 2°C - 4°C for adequate protection inside the container till implantation. Protection can safely extend to over 48 hours.

LIVER

RECOMMENDED SOLUTION TEMPERATURE: 0-4°C.

For a correct infusion of isolated organs, the following instructions should be followed:

INFUSION VOLUME: If liver, pancreas and kidneys are to be protected as a whole in the donor, an infusion volume of 150-200 ml of Servator® H solution/kg bodyweight will be necessary. In a patient of about 70-80 kg of bodyweight, this multi organ protection implies a infusion volume of 8-12 liters of Servator® H solution.

INFUSION PRESSURE: Gravity infusion is performed by placing the container 1 m above the heart level.

INFUSION TIME: Under these dosage and pressure conditions, infusion time should be 10-15 minutes and, in no case less than 8 minutes.

ADDITIONAL MEASURES: The blood of an organ donor should be completely heparinized before proceeding to infusion with Servator® H solution.

INFUSION TECHNIQUE: Servator® H solution is infused through the infrarenal aorta or the iliac artery of the organ donor through a dedicated infusion tube (air-free system). Simultaneously to the beginning of gravity

infusion, the surgeon opens the vena cava in the donor's abdomen. This allows the solution to flow out unhindered.

The entire volume of the solution is administered through the abdominal aorta, so that all abdominal organs are protected. The biliary passages - either inside or outside the body - should be thoroughly rinsed out with minimum 100 ml of cold Servator® H solution. Rinsing is usually performed using calibrated catheters.

If only the liver or part of the liver (e.g., in case of a living donor) is to be removed without any other organ, the volume flushed should be proportionally reduced. Infusion time should not be, under any circumstances, less than 8 minutes; generally it should be 10-15 minutes. In this case an adequate infusion both of the arterial and the portal vein circulation should be ensured.

For a correct preservation of organs for transplantation, the following instructions should be followed:

TRANSPLANTATION: After being removed, the liver should be immersed in cold Servator® H solution for transportation. The organ should be completely covered by the solution.

It is generally agreed that ischemic time should normally not exceed 12-15 hours. If the liver is to be operated ex situ (e.g., for enucleation of a tumor), it should be stored in cold solution during the entire procedure. After completion of the so-called "bench procedure", it should be immediately autotransplanted.

PANCREAS

For protection of the pancreas, the considerations reported in paragraph "Liver" apply.

This means that the indications concerning all critical parameters (solution temperature, volume, pressure and infusion time) given in that paragraph are valid also for the pancreas.

VEIN AND/OR ARTERY TRANSPLANT

The vein for transplant (usually part of the great saphenous vein) and/or the artery for transplant (usually part of the internal thoracic artery) are cooled and stored in 50- 100 ml of Servator® H solution at 5° - 8°C. The segments of vein or artery are implanted after removal from the solution.

MULTI ORGAN PROTECTION

The basic infusion techniques are now largely standardized and described in the appropriate surgery textbooks. Regarding the infusion technique, gravity infusion with a infusion system of the largest possible caliber is worldwide accepted. The Servator® H solution has an extraordinarily low viscosity even at temperatures close to 0°C.

Consequently large volumes of solution at low temperature and moderate pressure can be administered in accordance with the infusion technique.

Multi organ protection is not linked to a definite volume of solution but to a short and limited infusion time of 8 - 10 minutes. Administration of large volumes of Servator® H solution at low temperature (0° - 4°C) allows a rapid and efficient cooling and consequent protection of the organs within the required time.

TRANSPORT OF A DONOR ORGAN

The hypothermic storage protocol differs from one hospital to another, but the "triple bag technique" has now been widely adopted internationally. The organ removed from the donor is packed in a specially designed sterile bag containing ice-cold Servator® H solution. The organ must be completely covered by the solution. The bag is sealed with an adhesive strip or similar device and is then placed in a second container which is also filled with Servator® H solution, to avoid presence of trapped air that might affect insulation and cooling. The organ, adequately protected by this double packing, is placed in a sterile plastic container, sealed with a safety lid. The sealed container is then placed in the transport case containing ice to keep the proper temperature conditions during the transfer journey. Documents containing information referring to the donor, the laboratory analysis and some blood samples from the donor are also included. The transfer of the organ in the Servator® H solution should be the fastest possible.

SAFE STORAGE CONDITIONS: The solution intended for transport must be kept at 2°C - 4°C.

INCORRECT USE AND OVERDOSAGE: The entry of large volumes of solution into the systemic circulation may lead to an excessive circulating volume with consequent electrolytic imbalance (hypocalcemia, hyponatremia, hypermagnesemia, hyperkalemia). An adequate monitoring of serum electrolytes is recommended. The

complete inactivation of the myocardium renders it susceptible to overstretching. Care must therefore be taken to ensure adequate venting of the ventricle.

Do not exceed recommended infusion volumes and pressures.

Special precaution should be taken when dealing with infants' and children's hearts.

The use of the Servator® H solution at wrong temperatures may cause the so-called "calcium paradox" and lead to destruction of myocardial cells once the organ is implanted.

Consequently, the solution must always be cooled before use as indicated in the operative protocols.

ADVERS REACTIONS: None known at present.

WARNINGS AND PRECAUTIONS REGARDING THE PRODUCT VALIDITY: Verify the expiration date on the packaging. The expiration date refers to the product in its integral container. ATTENTION: do not use the solution after the expiration date.

Do not use the product, even before the expiration, if the solution has an intense yellow coloration equal to, or stronger than Y5 (see coloration scale of Eur. Ph. current edition).

Use the solution only if clear and without visible particles.

STORAGE CONDITIONS: Store in refrigerator (2°C - 8°C) and protect from light. The solution is sterile and is intended for one single and continuous administration.

Discard any residue to avoid risk of contamination due to loss of sterility.

MEDICAL DEVICE FORM AND CONTENT: PVC-free bag or glass bottle.

SPECIAL PRECAUTIONS FOR DISPOSAL OF UNUSED RESIDUES: Unused residues of the solution should be disposed of in conformity with the local rules in force.

GENERAL INFORMATION AND MECHANISM OF ACTION

Servator® H solution prolongs ischemic tolerance of the protected organ mainly through two mechanisms:

1. The electrolytic composition of Servator® H solution prevents the activation of energy consuming processes by reducing to the minimum levels the demand of energy by the organ.
2. Anaerobic energy production is limited by the increasing inhibition of glycolysis due to decrease of pH consequent to accumulation of lactic acid. The buffer histidine/histidine HCl delays the decrease of pH in the organ tissues during ischemia, thus increasing the proportion of anaerobic glycolytic energy produced.

Potassium hydrogen 2-oxopentandioate acts as substrate for production of aerobic energy.

Tryptophan is claimed to have a membrane protective action.

Mannitol is associated to prevention of cell swelling.

The overall osmolarity of the solution is slightly lower than the osmolarity of the plasma and the intracellular space.

TOXICOLOGICAL PROPERTIES

The inflow of large volumes of solution into the systemic circulation can lead to excessive circulating volume and consequent electrolyte unbalance (hyponatremia, hypocalcemia, hyperkalemia, hypermagnesiumemia).

Tryptophan and histidine may reach elevated plasma levels during the first 24 hours. Notwithstanding, no adverse effects on metabolism have been observed.

PHARMACOKINETICS

Depending on the nature and duration of the operation, the operative methods and the weight of the patient, the volume of solution entering the systemic circulation should range from 0.1 to 3.0 liters.

Potassium hydrogen 2-oxopentandioate is metabolized mainly through the citric acid cycle.

Histidine and tryptophan are metabolized mainly in the liver, though partly also excreted by the kidneys.

Mannitol is eliminated unchanged by the kidneys.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP): the document is available at the following URL:
<https://cloud.salfspa.it:8443/index.php/s/SSCP.ServatorH>

SERIOUS INCIDENT REPORTING: the user and/or patient are recommended to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

MANUFACTURER: S.A.L.F. S.p.A. LABORATORIO FARMACOLOGICO via Marconi, 2 - 24069 Cenate Sotto (BG)
Italy

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FABRICANTE: S.A.L.F. S.p.A. LABORATORIO FARMACOLOGICO via Marconi, 2 - 24069 Cenate Sotto (BG) Itália