# SERVATOR® M

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Data di revisione del testo: Marzo 2023

# SERVATOR® M

SOLUTION FOR PRESERVATION OF ORGANS

STERILE APYROGENIC

NON-INJECTABLE - DO NOT USE FOR INTRAVENOUS INFUSION

**REF** : HTKE0810GPO

## DEVICE TYPE:

CLASS III MEDICAL DEVICE.

## QUALITATIVE AND QUANTITATIVE COMPOSITION

COMPONENT	g/L	mM/L
Adenine (free base)	0.68	5
Calcium chloride dihydrate	0.068	0.5
Glucose monohydrate equal to anhydrous glucose	2 1.80	10
Reduced glutathione	0.92	3
HEPES free acid	2.38	10
Hydroxyethyl starch	50	NA
Magnesium gluconate hydrate equal to anhydrous Magnesium gluconate	1.18 1.13	5
Mannitol	5.4	30
Monobasic potassium phosphate	3.4	25
Ribose, D(-)	0.75	5
Sodium gluconate	17.45	80
Sodium hydroxide	0.70	NA
Water for injection preparations q.s. for 1000 ml		

## PHYSICAL PROPERTIES:

The solution has a theoretical osmolality of approximately 300 mOsmol/Kg, with a final concentration of sodium ions: 100 mEq/l and a contraction of potassium ions: 25 mEq/L.

pH: approximately 7.4 at room temperature.

## DESCRIPTION

SERVATOR® M is a sterile, apyrogenic, clear solution, transparent to pale yellow in colour, in PVC-free bag of 1000 ml. The bag is wrapped in an overbag. The bag can be stored at a temperature between +2°C and +25°C. It must be cooled to between +2°C and +8°C before use.

## INTENDED USE

SERVATOR® M is a sterile, pyrogen-free solution, indicated for rinsing and hypothermic storage of the kidney, isolated at the time of removal from the donor, prior to storage, transportation and possible transplantation into the recipient.

The medical device may only be used by suitably trained medical staff and according to established operating protocols.

#### **INSTRUCTIONS FOR USE:**

Before use, cool the Servator® M solution to a temperature between +2°C and +8°C. The ideal temperature is +5°C.

Remove the overbag.

The solution must not be filtered. Check that the solution is clear and colourless or pale yellow and free of visible particles. Check the tightness of the primary container. If a leak is detected, the bag cannot be used. The cooled solution does not come into direct contact with the donor patient but with the suitably isolated organ.

Servator M® must be flushed from the donor organ at the time of implantation.

#### **FLUSHING:**

Before transplantation in the recipient, the kidney must be rinsed so as to eliminate any residual quantity of Servator M solution using appropriate flushing solutions. Flushing is to be conducted in compliance with the transplantation center's current procedures. The flushing operations are to be implemented by competent medical personnel.

#### **STORAGE CONDITIONS**

Store in the original container at a temperature between +2°C and +25°C.

Avoid exposure to heat sources. Do not freeze. Protect from light sources. Cool the solution before use. Remove the overbag immediately before using the medical device. The ideal temperature for use of the medical device is +5°C; however, use between +2°C and +8°C is acceptable.

KEEP OUT OF REACH AND SIGHT OF CHILDREN.

#### **WARNINGS AND PRECAUTIONS BEFORE USE:**

DO NOT USE FOR DIRECT INJECTION OR INTRAVENOUS AND INTRA-ARTERIAL INFUSION.

DO NOT USE FOR RECONSTITUTION OF BLOOD VOLUME OR REPLACEMENT OF ELECTROLYTES OR AMINO ACIDS.

Systemic administration may lead to life-threatening side effects.

CAUTION: Failure to remove any residues of the preservation solution may result in hypersensitivity to one of its components. Hydroxyethyl starch has caused hypersensitivity reactions in some patients when used as an intravenous infusion directly. Medical staff must be prepared to respond to possible reactions to components or additives.

#### **ADVERSE REACTIONS**

No undesirable effects are known when the device is used as directed. The recipient may experience hypersensitivity to the components of the solution if the kidney flushing procedure prior to implantation is not carried out correctly.

#### **WARNINGS AND PRECAUTIONS IN CONNECTION WITH THE VALIDITY OF THE DEVICE**

See the expiry date specified on the package. Batch and expiry date are marked on each bag. The expiry date refers to a product in an intact package.

**CAUTION:** Do not use the device after the indicated expiry date.

Although within its validity (medical device not expired), the solution must not be used if any visible particles, precipitate or contamination are detected.

#### **SPECIAL PRECAUTIONS FOR DISPOSAL AND HANDLING OF THE DEVICE**

The unused medicinal product and the waste resulting from such a drug must be disposed of in compliance with the local regulations in force.

#### **NATURE AND CONTENTS OF THE PACKAGE**

1000 ml PVC FREE bag Cardboard box containing 10 bags.

**SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP):** The document is available at the following URL:

<https://cloud.salfspa.it:8443/index.php/s/SSCP.ServatorH>

**SERIOUS INCIDENT REPORTING:** The user and/or patient are recommended to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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S.A.L.F. S.p.A. LABORATOIRE PHARMACOLOGIQUE